



CORFLEX®

Designed for life.

CORFLEX.COM | 800 426 7353

Lower Leg & Ankle Pneumatic Walking Splint

Application Instructions

Fitting and Application:

1. Release straps from liner and remove from metal uprights. Place injured foot into liner and fit heel up against back of liner. Close forefoot cover over the foot. Secure the upper calf liner with hook and loop closure. Insert optional ankle pads now if desired.



2. Place the foot with liner into walking splint. Thread the strap closest to the toes through plastic buckle and secure onto itself. Repeat this process for all remaining straps working your way towards the knee. Center the uprights in line with the leg and remove plastic sleeves and engage the hook to body of the liner.



3. Test for desired comfort level, do not overtighten straps. To inflate splint liner turn the inflate/deflate knob clockwise and press and release blue air bulb several times until liner is providing comfortable support.

DO NOT OVERINFLATE.



CORFLEX®

Designed for life.

CORFLEX.COM | 800 426 7353

Lower Leg & Ankle Pneumatic Walking Splint

Application Instructions

Fitting and Application:

1. Release straps from liner and remove from metal uprights. Place injured foot into liner and fit heel snugly in posterior section. Close forefoot cover over the foot. Secure the upper calf liner with hook and loop closure. Insert optional malleolus pads now if desired.



2. Place the foot with liner into walking splint. Thread the strap closest to the toes through plastic buckle and secure onto itself. Repeat this process for all remaining straps working your way towards the knee. Center the uprights in line with the leg and remove plastic sleeves and engage the hook to body of the liner.



3. Test for desired comfort level, do not overtighten straps. To inflate splint liner turn the inflate/deflate knob clockwise and press and release blue air bulb several times until liner is providing comfortable support.

DO NOT OVERINFLATE.



APPLICATION INSTRUCTIONS (continued)

4. To release air from the splint liner, turn the inflate/deflate knob counterclockwise several times. When reinflating, turn knob clockwise until snug and begin pumping. If at any time the liner feels too tight while wearing product, immediately release air by turning inflate/deflate knob counterclockwise. Consult your medical professional immediately.



Care Instructions

The Corflex Pneumatic Walking Splint liner may be hand washed using warm water and detergent for delicate fabrics. Repeatedly rinse and line dry. Do not use bleach or place in a washer or dryer. Wipe plastic shell with damp cloth, you may use disinfectant or mild soap, and wipe dry.

Caution

DO NOT OVERINFLATE. Release air immediately if liner feels too tight. This product has been designed and tested based on single patient usage. Since it is a single use device it should never be used by multiple patients. If additional pain or symptoms occur while using this device, seek medical attention.

Indications

For post injury or post-operative application to grade 1, 2, or 3 ankle sprains, and stable or internally fixed fractures of the foot or ankle.

Contraindications

For unstable fractures or for fractures of the proximal tibia or fibula.

Warranty

Products manufactured by Corflex are warranted to the original purchaser against defects in workmanship and material for the life of the product. The foregoing warranty shall not apply to normal wear and tear; or if the product has been subjected to misuse or misapplication the warranty will be voided at the discretion of Corflex.

To Reorder

Pneumatic Walking Splints

	XS	S	M	L	XL
Lower Leg	79-1079	79-1080	79-1081	79-1082	79-1083
Ankle	79-1089	79-1090	79-1091	79-1092	79-1093

Pneumatic Walking Splints Replacement Liners

	XS	S	M	L	XL
Lower Leg	79-1179	79-1180	79-1181	79-1182	79-1183
Ankle	79-1189	79-1190	79-1191	79-1192	79-1193

Sizing Guide	Mens	Women
X-Small	4 to 4 ½	3 to 5
Small	5 to 6	5 ½ to 7
Medium	6 ½ to 10	7 ½ to 11
Large	10 ½ to 12	11 ½ to 13
X-Large	12 ½+	14+

Measured by Shoe Size.



APPLICATION INSTRUCTIONS (continued)

4. To release air from the splint liner, turn the inflate/deflate knob counterclockwise several times. When reinflating turn knob clockwise until snug and begin pumping. If at any time the liner feels too tight while wearing product, immediately release air by turning inflate/deflate knob counterclockwise. Consult your medical professional immediately.



Care Instructions

The Corflex Pneumatic Walking Splint liner may be hand washed using warm water and detergent for delicate fabrics. Repeatedly rinse and line dry. Do not use bleach or place in a washer or dryer. Wipe plastic shell with damp cloth, you may use disinfectant or mild soap, and wipe dry.

Caution

DO NOT OVERINFLATE. Release air immediately if liner feels too tight. This product has been designed and tested based on single patient usage. Since it is a single use device it should never be used by multiple patients. If additional pain or symptoms occur while using this device, seek medical attention.

Indications

For post injury or post-operative application to grade 1, 2, or 3 ankle sprains, and stable or internally fixed fractures of the foot or ankle.

Contraindications

For unstable fractures or for fractures of the proximal tibia or fibula.

Warranty

Products manufactured by Corflex are warranted to the original purchaser against defects in workmanship and material for the life of the product. The foregoing warranty shall not apply to normal wear and tear; or if the product has been subjected to misuse or misapplication the warranty will be voided at the discretion of Corflex.

To Reorder

Pneumatic Walking Splints

	XS	S	M	L	XL
Lower Leg	79-1079	79-1080	79-1081	79-1082	79-1083
Ankle	79-1089	79-1090	79-1091	79-1092	79-1093

Pneumatic Walking Splints Replacement Liners

	XS	S	M	L	XL
Lower Leg	79-1179	79-1180	79-1181	79-1182	79-1183
Ankle	79-1189	79-1190	79-1191	79-1192	79-1193

Sizing Guide	Mens	Women
X-Small	4 to 4 ½	3 to 5
Small	5 to 6	5 ½ to 7
Medium	6 ½ to 10	7 ½ to 11
Large	10 ½ to 12	11 ½ to 13
X-Large	12 ½+	14+

Measured by Shoe Size.



669 EAST INDUSTRIAL PARK DRIVE
MANCHESTER, NH 03109
CORFLEX.COM | TOLL FREE: 800 426 7353



669 EAST INDUSTRIAL PARK DRIVE
MANCHESTER, NH 03109
CORFLEX.COM | TOLL FREE: 800 426 7353