

# ReAlign Spinal Orthoses



**CORFLEX®**  
Designed for life.

CORFLEX.COM

**One Stocking Sku/Style - Universal Sizing – 30”- 50” & up to 75”  
with Extension Panels.**



Orthoses are designed to help stabilize and decompress the spine in patients recovering from surgery or injury. Rigid anterior and posterior panels limit flexion and extension and help protect and support the vertebrae and surrounding tissues. Hand loop and dual Power-Pull Closures™ make the brace easy to apply and deliver patient-controlled compression to the abdominal cavity to help unload the spine and support musculature.

NEW Cord Wrap Cord Winder™ allows precise length adjustment to properly position pull tabs across the universal size range. Step-down design allows for long-term patient value, simply remove back/side panels for a lightweight, streamlined design.

- Rigid anterior and posterior panels stabilize the spine and can be adjusted with heat to achieve a custom fit; hand loop and Power-Pull Closures allow patients with limited upper extremity strength and dexterity to independently apply the brace
- Indicated for lower back pain, sprains or strains, spinal stenosis, spondylosis, spondylolysis, spondylolisthesis, degenerative disc disease, bulging or herniated discs, burst fractures, facet syndrome and post-operative stabilization following laminectomy, foraminotomy, IDET procedures, discectomy, fusion or decompression
- Constructed of nylon, polyester/foam laminates, molded pull tabs, heat moldable ABS posterior panel, heat moldable HDPE anterior/lateral panels and hook and loop closures



Convenient cold/hot therapy pocket and pouch - fits 5 1/2" x 11" gel pack and attaches easily with hook and loop.



Easy to adjust Cord Winder allows for proper cord length especially at smaller sizes.



**Product Number**

**UNIV**

ReAlign Lumbar Orthosis (LO)

35-3430

Sizing 30" - 50"

Measurement taken from waist or hip circumference.

Suggested Billing Code L-0627, (Custom Fit), L-0642 (OTS)

Lo sag ri an/pos pnl pre cst, Lo sag ri an/pos pnl pre ots



**Product Number**

**UNIV**

ReAlign Lumbosacral Orthosis (LSO)

35-3440

Sizing 30" - 50"

Measurement taken from waist or hip circumference.

Suggested Billing Code L-0631 (Custom Fit), L-0648 (OTS)

Lso sag r an/pos pnl pre cst, Lso sag r an/pos pnl pre ots



**Product Number**

**UNIV**

ReAlign Lumbosacral Plus (LSO Plus)

35-3450

Sizing 30" - 50"

Measurement taken from waist or hip circumference.

Suggested Billing Code L-0637 (Custom Fit), L-0642 (OTS)

Lso sc r ant/pos pnl pre cst, Lo sag ri an/pos pnl pre ots



ReAlign Extension Kit

**Accessories**

**UNIV**

ReAlign Extension Kit, 16"

35-3912

Cryotherm Gel Pouch Only

72-7008

Cryotherm Pouch with Gel

72-7009

Gel Pack, 5 1/2" x 11"

72-7001

Gel Pack, 5 1/2" x 11", cs/10

72-7000



CORFLEX | 669 EAST INDUSTRIAL PARK DRIVE | MANCHESTER, NH 03109  
CORFLEX.COM | TOLL FREE: 800 426 7353 | FAX: 603 623 4111