



CORFLEX®

Designed for life.

NEW!

Contender™ Post-Op Knee Brace

The Contender Post-Op Knee Brace features a simple and precise adjustment system to minimize fit time after surgery, and at 27 oz., the brace is one of the lightest post-op knee braces available.



Contender™ Post-Op Knee Brace The brace that stays in place.

ZipTrak™ System

Adjust overall length as well as individual cuff positions in 1/8" increments

Trim to Fit Straps

Simply cut to length and apply alligator-style hook to ends

Drop Lock Slide Button

Six locking positions to quickly immobilize the brace from -10° to 40°

Drop Lock Degree Indicator

Easily visible by the patient during wear to promote protocol compliance

ZipTrak™ Lock

Allows quick, tool-free adjustment of brace length from 18"-26" and location of intermediate cuffs to avoid surgical sites

Padded Liner

Slip resistant foam liners help maintain proper brace position

Ankle Cuff

Smaller in size and highly flexible; sits just proximal to malleolus to further prevent migration

Length Indicators

Graduated markings on the rails facilitate quick fitting following surgery

Metal Uprights

Custom bendable, aircraft aluminum uprights are notched to assist in contouring

FlexLok™ & ROM Buttons

Finger touch controls adjust range of motion from -10° to 110° and feature a unique FlexLok band to help prevent patient tampering

QuikZip™ Buckles

Allow for easy application and micro adjustment of strap for a comfortable and secure fit

X-Flex™ Cuff System

Flexes to hug the leg, significantly reducing brace migration

Contender™

The brace that stays in place.

Available in two styles



Lite Version



Full Foam Version

- > **27 oz.** Notably lighter than competitive braces without sacrificing strength
- > **ZipTrak™ System** Facilitates length adjustment from 18"-26" in just seconds for a precise fit
- > **X-Flex™ Cuff System** Conforms to hug the leg, significantly reducing brace migration
- > **QuikZip™ Buckles** Allow for easy application and micro adjustment for a comfortable and secure fit

Find out more at
corflex.com/postopknee

Indicated for adjustable range of motion control following ACL, PCL, MCL and LCL surgeries, meniscal repairs, patella realignment, regenerative chondroplasty, stable femoral fractures, total knee replacements and high tibial osteotomies

Product Number

	Brace	Replacement Liner
Lite	75-7500	75-7501
Full Foam	75-7550	75-7551

Sizing Guide

Fits up to 34" thigh circumference

Length Adjustable 18"-26"

*Measurement taken from leg circumference
6" above mid patella*

Suggested Billing Codes

Knee Orthosis, Adjustable Knee Joints
L-1832 (Custom-Fit)
or L-1833 (OTS)



CORFLEX®

Designed for life.

800 426 7353

CORFLEX.COM

669 EAST INDUSTRIAL PARK DRIVE / MANCHESTER, NH 03109

FAX 603 623 4111 / EMAIL SERVICE@CORFLEX.COM



ALL PRODUCTS MANUFACTURED BY
CORFLEX ARE LATEX FREE.