

Cryotherm™ and Cryo Pneumatic Rehabilitation Products

A complete line of anatomically contoured wraps designed to provide compression and concentrated cold/hot therapy for various parts of the body. Both styles are ideal for post-surgical rehab, inflammation, physical therapy, occupational therapy and athletic training.



CORFLEX®

Designed for life.

CORFLEXGLOBAL.COM




Cryo Pneumatic Knee Splint

- Anatomical knee wrap designed to combine the benefits of pneumatic compression with cold therapy; features pump bulb that detaches once desired compression is achieved
- Indicated for post-op application or as part of a physical therapy or athletic training program; helps diminish swelling/edema, decrease pain and reduce muscle spasms
- Manufactured from plush foam laminate with contact closure straps, removable gel pack with tricot lining and non-toxic gel; additional universal extension straps available

Product Number	0 Gels	1 Gel	2 Gels
	72-7170	72-7171	72-7172

Sizing Guide	Univ
	Fits up to 24" circumference Measurement taken from leg circumference 4" above mid-patella

Accessories	
Gel Pack, each	72-7179
Gel Pack, 10/bx	72-7179-010
14" Universal Ext Straps (10 Pack)	72-7190
Pump & Hose Set	72-7165



Gel Pack
Removable, contoured



Valve
On/off valve allows for easy inflation and deflation



Bulb
Features button for quick deflation





Cryo Pneumatic Knee Orthosis w/Hinge

- Anatomical knee wrap that features bi-lateral hinges providing medial/lateral support; movable hinges may be contoured and feature a built in stop at full extension; includes two extension straps to facilitate hinge placement
- Indicated for post knee arthroscopy, meniscal tears, osteoarthritis pain and swelling and patellofemoral pain
- Manufactured from plush foam laminate with hook and loop closures; reusable gel packs feature soft fabric lining and non-toxic gel; includes two aluminum hinges and extension straps

Product Number	0 Gels	1 Gel	2 Gels
	72-7270	72-7271	72-7272

Sizing Guide	Univ
	Fits up to 24" circumference Measurement taken from leg circumference 4" above mid-patella

Accessories	
Gel Pack, each	72-7179
Gel Pack, 10/bx	72-7179-010
14" Universal Ext Straps (10 Pack)	72-7190
Pump & Hose Set	72-7165

Suggested Billing Code L-1810 (Custom Fit), L-1812 (OTS)



Cryo Pneumatic Knee Orthosis w/R.O.M. Hinge

- Features ROM Hinges with malleable aluminum uprights offering full flexion and extension control with settings at 0°, 15°, 30°, 60° and 90°
- Indicated following ACL repair, post total joint surgery, meniscal tears, patella realignment or during rehabilitation
- Manufactured from plush foam laminate with hook and loop closures; reusable gel packs feature soft fabric lining and non-toxic gel; includes two aluminum ROM hinges with plastic covers and extension strap

Product Number	0 Gels	1 Gel	2 Gels
	72-7370	72-7371	72-7372

Sizing Guide	Univ
	Fits up to 24" circumference Measurement taken from leg circumference 4" above mid-patella

Accessories	
Gel Pack, each	72-7179
Gel Pack, 10/bx	72-7179-010
14" Universal Ext Straps (10 Pack)	72-7190
Pump & Hose Set	72-7165

Suggested Billing Code L-1832 (Custom Fit), L-1833 (OTS)





Cryotherm™ Knee Wrap

- Economical wrap designed to provide compression and concentrated cold/hot therapy for the knee; available with 2 pockets or 3 pockets for circumferential coverage
- Indicated for post-op use, arthritis, strains or sprains and for integration into rehab program
- Manufactured from soft breathable foam laminate with contact closure and pockets to hold reusable gel packs

Product Number	0 Gels	2 Gels	3 Gels	4 Gels	6 Gels
2 Pockets	72-7020	72-7022		72-7024	
3 Pockets	72-7080		72-7083		72-7086

Sizing Guide	Univ
	Fits up to 24" circumference Measurement taken from leg circumference 4" above mid-patella

Accessories	
Gel Pack, 5 1/2" x 11", each	72-7001
Gel Pack, 5 1/2" x 11", 10/bx	72-7000



Cryotherm™ Back Wrap

- Economical wrap designed to provide compression and concentrated cold/hot therapy for the back
- Indicated for post-op use, arthritis, sprains/strains and for integration into rehab program
- Manufactured from a soft breathable foam laminate with double-pull contact closure and pockets to accommodate reusable gel packs

Product Number	0 Gels	2 Gels	4 Gels
	72-7040	72-7042	72-7044

Sizing Guide	Univ
	Fits up to 48" chest circumference

Accessories	
Gel Pack, 5 1/2" x 11", each	72-7001
Gel Pack, 5 1/2" x 11", 10/bx	72-7000



Cryo Pneumatic Ankle Wrap

- Anatomical ankle wrap designed to combine the benefits of pneumatic compression with cold therapy; features pump bulb that detaches once desired compression is achieved
- Indicated for post-op application, severe ankle strains/sprains, and as part of a physical therapy or athletic training program; helps diminish swelling/edema, decrease pain and reduce muscle spasms
- Manufactured from plush foam laminate that is hook compatible; removable gel pack with tricot lining; non-toxic gel

Product Number	0 Gels	1 Gel	2 Gels
	72-7210	72-7211	72-7212
Sizing Guide	Univ		
	Fits up to 16" ankle circumference Measurement is circumference from base of heel around ankle		
Accessories			
Gel Pack, each		72-7219	
Gel Pack, 10/bx		72-7219-010	



Cryotherm™ Ankle Wrap

- Economical wrap designed to provide compression and concentrated cold/hot therapy for the ankle
- Indicated for post-op use, arthritis, strains or sprains and as part of rehab program
- Manufactured from a soft breathable foam laminate with contact closure and pockets to hold reusable gel packs

Product Number	0 Gels	2 Gels	4 Gels
	72-7140	72-7142	72-7144
Sizing Guide	Univ		
	Fits up to 14" ankle circumference Measurement is circumference from base of heel around ankle		
Accessories			
Gel Pack, 3 1/2"x 11", each		72-7003	
Gel Pack, 3 1/2"x 11", 10/bx		72-7002	



Cryo Pneumatic Shoulder

- Anatomical shoulder wrap designed to combine the benefits of pneumatic compression with cold therapy; features pump bulb that detaches once desired compression is achieved
- Indicated for post-op application, AC joint separations, bursitis, sports injuries and as part of a physical therapy or athletic training program; helps diminish swelling/edema, decrease pain and reduce muscle spasms
- Manufactured from plush foam laminate with contact closure straps, removable gel pack with tricot lining and non-toxic gel; additional universal extension straps available

Product Number	0 Gels	1 Gel	2 Gels
	72-7180	72-7181	72-7182

Sizing Guide	Univ
	Fits up to 48" chest circumference

Accessories	
Gel Pack, each	72-7189
Gel Pack, 10/bx	72-7189-010
14" Universal Ext Straps (10 Pack)	72-7190
Pump & Hose Set	72-7165

Suggested Billing Code L-3670



Cryotherm™ Shoulder Wrap

- Designed to contour to the shoulder with pockets that hold gel packs to provide benefits of cold/hot therapy
- Indicated for post-op use, treatment of arthritis or as part of a physical therapy or rehab program
- Manufactured from soft breathable foam laminate with contact closures and pockets to hold reusable gel packs

Product Number	0 Gels	2 Gels	4 Gels
	72-7030	72-7032	72-7034

Sizing Guide	Univ
	Fits up to 48" chest circumference

Accessories	
Gel Pack, 5 1/2" x 11", each	72-7001
Gel Pack, 5 1/2" x 11", 10/bx	72-7000





Cryo Pneumatic Elbow Wrap

- Anatomical elbow wrap designed to combine the benefits of pneumatic compression with cold therapy; features pump bulb that detaches once desired compression is achieved
- Indicated for post-op application, epicondylitis, cubital tunnel syndrome and as part of a physical therapy or athletic training program; helps diminish swelling/edema, decrease pain and reduce muscle spasms
- Manufactured from plush foam laminate that is hook compatible; removable gel pack with tricot lining; non-toxic gel

Product Number	0 Gels	1 Gel	2 Gels
	72-7210	72-7211	72-7212
Sizing Guide	Univ		
	Fits up to 16" elbow circumference		
Accessories			
Gel Pack, each		72-7219	
Gel Pack, 10/bx		72-7219-010	



Cryotherm™ Elbow Wrap

- Economical wrap designed to provide compression and concentrated cold/hot therapy for the elbow
- Indicated for post surgical application, arthritis, strains or sprains and for integration into rehab program
- Manufactured from soft breathable foam laminate with contact closure and pocket to hold reusable gel packs

Product Number	0 Gels	1 Gel	2 Gels
	72-7050	72-7051	72-7052
Sizing Guide	Univ		
	Fits up to 18"bicep circumference Measurement taken at widest part of bicep		
Accessories			
Gel Pack, 5 1/2"x 11", each		72-7001	
Gel Pack, 5 1/2"x 11", 10/bx		72-7000	



Cryo Pneumatic Wrist Support

- Anatomical shoulder wrap designed to combine the benefits of pneumatic compression with cold therapy; features pump bulb that detaches once desired compression is achieved
- Indicated for post-op application, carpal tunnel syndrome, severe sprains or strains and as part of a physical therapy or occupational therapy program; helps diminish swelling/edema
- Manufactured from plush foam laminate that is hook compatible; removable gel pack with tricote lining; non-toxic gel; movable palmar stay

Product Number	0 Gels	1 Gel	2 Gels
Right	72-7240	72-7241	72-7242
Left	72-7250	72-7251	72-7252

Sizing Guide	Univ
	Fits up to 9" wrist circumference

Accessories	
Gel Pack, each	72-7259
Gel Pack, 10/bx	72-7259-010



Cryotherm™ Wrist Wrap

- Economical wrap designed to provide compression and concentrated cold/hot therapy for the wrist
- Indicated for post-op use, arthritis, wrist strains or sprains and for integration into rehab program
- Manufactured from a soft breathable foam laminate with contact closures and pockets to hold reusable gel packs

Product Number	0 Gels	1 Gel	2 Gels
	72-7150	72-7151	72-7152

Sizing Guide	Univ
	Fits up to 9" wrist circumference

Accessories	
Gel Pack, 3 1/2" x 11", each	72-7003
Gel Pack, 3 1/2" x 11", 10/bx	72-7002



Cryotherm™ Hip Wrap

- Designed to contour to the hip and provide mild compression with the benefits of cold/hot therapy; wrap can be positioned to provide coverage to the anterior, lateral or posterior hip anatomy
- Indicated post-surgically following Total Hip Replacement, Resurfacing, Revision or Arthroscopy; non-surgically following injections to the hip joint, arthritis, bursitis, tendonitis, muscle or tendon strains, labral tears, femoroacetabular impingement, piriformis syndrome, sciatica and hip dislocation/subluxation
- Manufactured from a soft breathable foam laminate with contact closures and pockets to hold two reusable gel pack

Product Number	0 Gels	2 Gels	4 Gels
	72-7060	72-7062	72-7064

Sizing Guide	Univ
	Fits up to 54" waist circumference and 36" thigh circumference

Accessories	
Gel Pack, 5 ½" x 11", each	72-7001
Gel Pack, 5 ½" x 11", 10/bx	72-7000



CORFLEX GLOBAL | 669 EAST INDUSTRIAL PARK DRIVE | MANCHESTER, NH 03109
CORFLEXGLOBAL.COM | TOLL FREE: 800 426 7353 | FAX: 603 623 4111