

CoolTex™ Ag 6" Wrist Splint

Low profile design and flexible support
ensure maximum patient comfort



CORFLEX®


Designed for life.

CORFLEXGLOBAL.COM



Our 6" CoolTex Ag Wrist Splint is designed to provide secure and comfortable immobilization of the wrist while preserving finger mobility. Single piece design minimizes seams to ensure maximum patient comfort while providing uniform compression. Wide elastic wrist strap provides additional compression that is easily controlled with a hook and loop closure. Splint features heat-moldable, plastic palmar stays for optimal stabilization. Elastic wrist strap provides additional support and compression.

Lightweight and breathable CoolTex Ag material allows for greater comfort during extended wear, helping wick moisture away from the skin. Material also contains silver yarns that provide antimicrobial protection to the fabric to help reduce odors and protect fibers from premature breakdown.



ADDITIONAL FEATURES:

- Ideal for all-day wear, support levels can be easily adjusted by removing stays or changing compression support at the wrist
- Indicated for Carpal Tunnel Syndrome, tendinitis, sprains or strains, soft tissue injury, wrist contracture and post-cast removal
- Constructed of breathable, moisture-wicking, antimicrobial CoolTex Ag liner; high performance brushed nylon laminate exterior; lightweight, moldable palmar stays

CoolTex™ Ag 6" Wrist Splint



Dual dorsal stays made from semi-rigid LDPE and HDPE plastic are removable to allow varied support based on activity.

Product Number Right	Product Number Left	Size	Measurement (Wrist Circumference)
74-5620	74-5630	XS	5" - 5 3/4"
74-5621	74-5631	S	5 3/4" - 6 1/2"
74-5622	74-5632	M	6 1/2" - 7 1/4"
74-5623	74-5633	L	7 1/4" - 8"
74-5624	74-5634	XL	8" - 8 3/4"
74-5625	74-5635	2XL	8 3/4" - 9 1/2"

Suggested Billing Code

L-3908 Wrist cock-up, non-molded



Lightweight, moldable palmar stays are contoured and removable for individual comfort level.

ALSO AVAILABLE

CoolTex™ Ag 6" Thumb Spica

The shorter length of our 6" CoolTex Ag Thumb Spica is ideal for use during daily activities. Removable stays allow for customized support of the thumb and wrist independently.

**ASK ABOUT OUR FULL LINE OF
COOLTEX AG PRODUCTS**



CORFLEX GLOBAL | 669 EAST INDUSTRIAL PARK DRIVE | MANCHESTER, NH 03109
CORFLEXGLOBAL.COM | TOLL FREE: 800 426 7353 | FAX: 603 623 4111