

CoolTex™ Ag 8" Thumb Spica with Strap

Silver-based antimicrobial lining provides continuous protection against odor causing microbes



CORFLEX®

Designed for life.

CORFLEXGLOBAL.COM



The CoolTex Ag Thumb Spica is designed to provide secure and comfortable immobilization of the wrist and thumb while preserving finger mobility. Single piece design with ulnar closure minimizes seams to ensure maximum patient comfort. Splint features aluminum removable and malleable palmar, thumb and dorsal stays for optimal stabilization. Elastic wrist strap provides additional support and compression.

The brace is manufactured from CoolTex Ag fabric to allow for easy flow of air, transfer of moisture and antimicrobial protection. Lightweight and breathable material allows for greater comfort during extended wear. Available in seven sizes, XS through 3XL in left and right.

Cool Airflow In **Heat and Perspiration Out**

Antimicrobial Silver Lining

CoolTex™ Ag
Freshness is in the Silver Lining

- A fighting layer of defense that continuously inhibits the growth of microbes, helping fabrics last longer and limiting odor formation.
- Breathable, open-mesh design promotes the transfer of heat and moisture away from the skin keeping you cool and dry.
- Nylon and spandex yarns combine to deliver superior compression and support while being extremely lightweight.

ADDITIONAL FEATURES:

- 8" splint provides compression and stabilization of the wrist and thumb in a unique, single-piece design for ultimate comfort
- Indicated for deQuervain's Syndrome, tendinitis, thumb strains or sprains and post-cast removal
- Manufactured from lightweight, breathable, antimicrobial CoolTex Ag material; elastic; aluminum stays

Suggested Billing Codes

- L-3807 (Custom Fit) WHFO, no joint, prefabricated
- L-3809 (OTS) WHFO without joints prefabricated

Product Number Right	Product Number Left	Size	Measurement (Wrist Circumference)
74-8600	74-8610	XS	4"-5"
74-8601	74-8611	S	5"-6"
74-8602	74-8612	M	6"-7"
74-8603	74-8613	L	7"-8"
74-8604	74-8614	XL	8"-9"
74-8605	74-8615	2XL	9"-10½"
74-8606	74-8616	3XL	10½"-12"



CORFLEX GLOBAL | 669 EAST INDUSTRIAL PARK DRIVE | MANCHESTER, NH 03109
CORFLEXGLOBAL.COM | TOLL FREE: 800 426 7353 | FAX: 603 623 4111