

# CoolTex™ Ag 13" Knee Sleeve with Stays

Open patella knee sleeve featuring CoolTex Ag fabric for breathable comfort



**CORFLEX®**

Designed for life.

CORFLEXGLOBAL.COM



The CoolTex Ag 13" Knee Sleeve with Stays features a slip-on design that provides support and compression to the knee and surrounding musculature. Brace features dual bilateral flex stays to help prevent sleeve rollover during activity.

The sleeve is manufactured from CoolTex Ag fabric to allow for easy flow of air, transfer of moisture and antimicrobial protection. Lightweight and breathable material allows for greater comfort during extended wear. Available in both open and closed popliteal, in eight sizes, XS through 4XL.

**Cool Airflow In**      **Heat and Perspiration Out**

**Antimicrobial Silver Lining**

**CoolTex™ Ag**  
Freshness is in the Silver Lining

- A fighting layer of defense that continuously inhibits the growth of microbes, helping fabrics last longer and limiting odor formation.
- Breathable, open-mesh design promotes the transfer of heat and moisture away from the skin keeping you cool and dry.
- Nylon and spandex yarns combine to deliver superior compression and support while being extremely lightweight.

#### ADDITIONAL FEATURES:


- 13" sleeve provides support and compression for the knee with flex stays to prevent roll-over; open patella aids in tracking
- Indicated for moderate knee strains or sprains, knee injuries, arthritis and post-op rehab
- Manufactured from antimicrobial CoolTex Ag, a lightweight, breathable neoprene alternative; plastic flex stays




# CoolTex™ Ag 13" Knee Sleeve with Stays



Open popliteal design promotes greater comfort, especially during activities that involve frequent knee flexion (also available in closed popliteal).



Brace tabs assist with application and adjustment during wear.



Dual flex stays are lightweight and help prevent sleeve rollover during activity.

Product Number		Size	Measurement		
Closed Popliteal	Open Popliteal		6" Above Mid-Patella	Mid-Patella	6" Below Mid-Patella
85-0152	85-0162	XS	16"-18"	12"-13"	11"-12"
85-0153	85-0163	S	17"-19"	13"-14"	12"-13"
85-0154	85-0164	M	18"-20"	14"-15"	13"-14"
85-0155	85-0165	L	19"-21"	15"-16"	14"-15"
85-0156	85-0166	XL	20"-22"	16"-18"	15"-16"
85-0157	85-0167	2XL	21"-23"	18"-20"	16"-17"
85-0158	85-0168	3XL	23"-26"	20"-22"	17"-18"
85-0159	85-0169	4XL	26"-29"	22"-24"	18"-19"



CORFLEX GLOBAL | 669 EAST INDUSTRIAL PARK DRIVE | MANCHESTER, NH 03109  
CORFLEXGLOBAL.COM | TOLL FREE: 800 426 7353 | FAX: 603 623 4111